

Specifications

Dokument Nr. / Document Nr. :	SOP AR-Ü08.4
Erste Veröffentlichungsdatum / First publication date:	25.09.2009
Revisionsdatum / Revision date:	11.05.2017
Revision Nr. / Revision Nr.:	03
Seite / Page:	1 / 1

Specifications

PRODUKT NAME: ACTODERM
PRODUCT NAME

ARTIKEL NR.: 06.2501.
ARTICLE NR

Prüfparameter Test Parameters	Test methode Test Method	Limit Limits
Aussehen Appearance	Ph. Eur (Visual inspection)	Hellblau Flüssigkeit Light Blue Liquid
Geruch Odor	Ph. Eur (Sensorial inspection)	Charakteristisch Characteristic
pH (Konz.) pH (Conc.)	Ph. Eur 2.2.3 01/2008:20203 (pH-meter)	Min: 6,0 Max: 8,0
Dichte (20°C) Density (20°C)	Ph. Eur 2.2.5 01/2008:20205 (Pycnometer)	Min: 0,854 g/ml Max: 0,895 g/ml

Menge des Wirkstoffes In 100 g des Produkts (EN: Amount of active substance in 100 g of the product)

2-Propanol (CAS No: 67-63-0)	Ph. Eur 2.2.28 01/2008:20228 GC	Min: 56,52 g % v/v 63	Max: 69,08 g % v/v 77
---------------------------------	------------------------------------	--------------------------	--------------------------

Qualitätskontrolle

Quality Control



Qualitätssicherung

Quality Assurance



Forschung & Entwicklung

Research & Development

